

## **How GoCheck Kids Works**

GoCheck Kids is an FDA-cleared, <u>Class Ile</u> medical device and meets the definition of an instrument-based ocular screener.

GoCheck Kids is a cloud-based system, where **image analysis** is done both **on the screening device (on-site) and remotely** on the GoCheck Kids' HIPAAsecure cloud:

GoCheck Kids Image Analysis

- When a photoscreening is taken, the image is analyzed on device (onsite) and analyzed on the remote GoCheck Kids mobile servers (GoCheck Kids cloud).
- After the image is analyzed both on the device and on remote GoCheck Kids mobile server, results are presented to the user on the device
- The image and results are immediately uploaded to our remote web servers and storage. Results and images can be analyzed remotely by GoCheck Kids on this web server.

If you have any questions, please contact us at accountmgmt@gobiquity.com.

## Mobile Application (on device image analysis) Remote Mobile Application Server (image analysis) Remote Web Application Server (image review) Database Server Encrypted MongoDB GoCheck Kids Image Quality Analysis